


DyeVert Clinical Update

December 2020




**be kind to
KIDNEYS**

DyeVert – Commitment to Clinical Outcomes

- Validation of DyeVert System's impact on outcomes
 - Over 15 clinical initiatives
 - 22 hospitals with over 70 users, and > 1,500 patients
- Stage 1: Performance outcomes in reducing AKI risk factors
 - Total contrast volume to patients
 - Total contrast relative to the patients baseline renal function
- Stage 2: Therapeutic outcomes in preventing AKI
 - Reduction of contrast delivered to patients
 - Implementing Threshold management

DYEVERT™
Contrast Reduction System **PLUS EZ**



DyeVert AKI Reduction Outcomes Data

Population Health Studies

Published and/or presented at Scientific Conferences



Control Group Studies

Three published at Scientific Conferences
One manuscript published in clinical journal




Economic Viability



Population Healthcare Quality Improvement


- 4 hospitals, 2016 – 2020
 - Quality improvement (QI) programs
 - Tracked outcomes through CathPCI AKI metric
- Patient centered program defined by clinical guidelines
 - Screening for risk
 - Volume management per SOC
 - Contrast monitoring and minimization
 - DyeVert Use in CKD/high risk patients
- **55% Reduction in AKI**
 - *35% - 84% AKI reduction range*



Author/Year	Baseline AKI Rate	Ending AKI Rate	Absolute AKI Reduction	Relative AKI Reduction
Castro et al, 2018	6.02%	3.90%	2.1%	35%
Thukral, 2019	11.69%	6.38%	5.3%	46%
Tucker et al, 2020	12.50%	2.04%	10.5%	84%
Cameron et al, 2020	11.45%	4.98%	6.5%	57%

DyeVert vs Control Group Comparative Data

- 4 hospitals, 2017 – 2020
- Clinical practice guidelines based program
 - Screening for risk
 - Volume management per SoC
 - Contrast monitoring and minimization
 - DyeVert Use in CKD/high risk patients
- Retrospective data abstraction of DyeVert & Control cases during the same time periods
- **51% Reduction in AKI**
 - 25% - 63% AKI reduction range



Author/Year	Control Group AKI Rate	DyeVert Group AKI Rate	Absolute AKI Reduction	Relative AKI Reduction
Sattar et al, 2018	16.20%	12.20%	4.0%	25%
Bunney et al, 2019	9.35%	3.45%	5.9%	63%
Kutschman, 2019	21.80%	9.40%	12.4%	57%
Briguori et al, 2020	19.00%	8.00%	11%	58%

Impact of a contrast media volume control device on acute kidney injury rate in patients with acute coronary syndrome

A Closer Look at Briguori, et al Propensity Matched Study comparing DyeVert vs. Control Group

- Purpose
 - Impact of DyeVert on contrast media in ACS patients
 - Impact of DyeVert on AKI rate in ACS patient
- Single-center, retrospective, two-arm study n=180
- Patients with Acute Coronary Syndrome
 - STEMI or high risk NSTEMI
 - Excluded ESRF w/Dialysis, recent CM exposure, and referrals for stroke center
- Hydration
 - IV 0.9% Sodium Chloride – all patients, except pulmonary congestion or CS
 - POSEIDON protocol intra procedure
- Propensity Matched Control
 - Age, Gender, LVEF, ACS type, BP, SCr, eGFR, Diabetes, PAD
- AKI analysis
 - SCr at baseline, 24, 48, and 72 hours
 - Calculated eGFR
 - CM Threshold Max = 3x eGFR

Impact of a contrast media volume control device on acute kidney injury rate in patients with acute coronary syndrome

Patients and Procedural Info:

- 90 patients included per group – control and DyeVert
- Primarily STEMI's – Control (81%) DyeVert (85%)
- Median Hydration volume – Control (1680 ml's) and DyeVert (1690 ml's)
 - 91% of controls hydrated > 960 ml's and 95.5% of DyeVert hydrated > 960 ml's
- ≥ 97% patients treated via radial approach
- All patients treated with PCI except for 3
- Mean CM volume – Control (130 ml's) and DyeVert (99 ml's)*
- CM Volume > 3x eGFR – Control (60%) and DyeVert (47%)
- DyeVert was never turned off under any circumstance due to inadequate image or device related reason

* $p < .001$

DyeVert vs Control Group

Outcomes

- AKI Rate in Control vs. DyeVert* Group
 - Control – 19% (17/90)
 - DyeVert – 8% (7/90)
- AKI reduction in DyeVert Group
 - 11% absolute vs. Control
 - 58% relative vs. Control
- Contrast reduction in DyeVert Group
 - 38% CMV savings vs. Control
- Reduced LOS in DyeVert** vs. Control Group

* $p < 0.047$ ** $p < 0.003$

Received 27 April 2020 | Revised 6 June 2020 | Accepted 27 June 2020
DOI: 10.1002/ccd.29136

ORIGINAL STUDIES

WILEY

Impact of a contrast media volume control device on acute kidney injury rate in patients with acute coronary syndrome

Carlo Briguori MD, PhD¹ | Marco Golino MD¹ | Nicola Porchetta MD¹ |
Mario Scarpelli MD¹ | Francesca De Micco MD¹ | Carmine Rubino RN² |
Amelia Focaccio MD¹ | Giuseppe Signoriello PhD²

¹Interventional Cardiology Unit, Mediterranean Cardiocentro, Naples, Italy
²Department of Mental Health and Preventive Medicine, Second University of Naples, Naples, Italy

Correspondence
Carlo Briguori MD, PhD, Interventional Cardiology, Mediterranean Cardiocentro, Via Orsini, 2, I-80121, Naples, Italy
Email: carlobriguori@ccicmediterranea.it

Abstract
BACKGROUND: The DyeVert™ system (Dyspey Medical Inc., Minnesota, MN) may reduce contrast media (CM) volume during coronary procedures while maintaining fluoroscopic image quality. Here, we assessed whether the use of the DyeVert system reduces acute kidney injury (AKI) rate in patients with acute coronary syndrome (ACS) undergoing invasive coronary procedures.
METHODS: ACS patients scheduled for coronary procedure from January 2017 to December 2019 were included. Two groups were identified: (a) Control group (n = 339), including patients in which a conventional manual injection syringe was used; and (b) DyeVert group (n = 112), in which CM injection was handled by the DyeVert™ system. A propensity score matching was performed to reduce the effect of treatment selection bias and potential confounders. In all cases, a low-osmolar, nonionic CM was administered. The primary objective was the rate of AKI, defined as a serum creatinine increase ≥ 0.3 mg/dL within 72 hr after CM exposure.
RESULTS: CM volume was higher in the Control group than in the DyeVert group (130 [120–188] mL vs. 99 [69–136] mL; $p < .001$). In the DyeVert group the mean percent CM volume saved was $38 \pm 13\%$. AKI occurred in 7/90 patients (8%) in the DyeVert group and in 17/90 (19%) patients in the Control group (odds ratio = 0.37; 95% confidence interval 0.14–0.95; $p = .047$).
CONCLUSIONS: This preliminary result suggests that CM volume reduction obtained by the DyeVert™ system is an effective strategy to prevent AKI in ACS patients undergoing invasive procedure.
KEYWORDS
acute coronary syndrome, acute kidney injury, contrast media

1 | INTRODUCTION
Acute kidney injury (AKI) is a common complication in patients suffering from acute coronary syndrome (ACS).^{1,2} This complication has been associated with higher early and late adverse events.^{1–4} The pathogenesis of AKI in ACS patients is multifactorial.^{1,4} Although some data suggest that it is unlikely that contrast media (CM) is a primary pathogenic factor responsible for renal dysfunction in ACS patient treated by percutaneous coronary intervention (PCI), avoiding excess CM volume has been advocated.^{4–6}
The DyeVert™ system (Dyspey Medical Inc., Minnesota, MN) is a novel device designed to reduce CM volume during coronary procedures, while maintaining fluoroscopic image quality.^{7,10} The aim of the present study is to assess whether the use of the DyeVert system is

Catheter Cardiovasc Interv. 2020;1–9.
wileyonlinelibrary.com/journal/ccd
© 2020 Wiley Periodicals LLC. | 1

DyeVert vs Control Group

Discussion:

- Patients with ACS at high risk for AKI
 - Ranges from 15 – 30%
- Majority of patients Stage 1 CKD at baseline
- Significant AKI Reduction in DyeVert group even with a high hydration regimen

Conclusion:

- *38% CM volume reduction with the DyeVert System is associated with significantly lower AKI rates vs. Control Group*

Received 27 April 2020 | Revised 6 June 2020 | Accepted 27 June 2020
DOI: 10.1002/ccd.20136

ORIGINAL STUDIES

Impact of a contrast media volume control device on acute kidney injury rate in patients with acute coronary syndrome

Carlo Briguori MD, PhD¹ | Marco Golino MD¹ | Nicola Porchetta MD¹ | Mario Scarpelli MD¹ | Francesca De Micco MD¹ | Carmine Rubino RN¹ | Amelia Focaccio MD¹ | Giuseppe Signorile PhD²

¹Interventional Cardiology Unit, Mediterraneo Cardiocenters, Naples, Italy
²Department of Mental Health and Preventive Medicine, Second University of Naples, Naples, Italy

Correspondence
Carlo Briguori MD, PhD, Interventional Cardiology, Mediterraneo Cardiocenters, Via Oratio, 2-80121, Naples, Italy.
Email: carlobriguori@clinicamediterranea.it

WILEY

Check for updates

Abstract

BACKGROUND: The DyeVert™ system (Isprey Medical Inc., Minnesota, MN) may reduce contrast media (CM) volume during coronary procedures while maintaining fluoroscopic image quality. Here, we assessed whether the use of the DyeVert system reduces acute kidney injury (AKI) rate in patients with acute coronary syndrome (ACS) undergoing invasive coronary procedures.

METHODS: ACS patients scheduled for coronary procedure from January 2017 to December 2019 were included. Two groups were identified: (a) Control group (n = 339), including patients in which a conventional manual injection syringe was used; and (b) DyeVert group (n = 112), in which CM injection was handled by the DyeVert™ system. A propensity score matching was performed to reduce the effect of treatment selection bias and potential confounders. In all cases, a low-osmolar, nonionic CM was administered. The primary objective was the rate of AKI, defined as a serum creatinine increase ≥ 0.3 mg/dL within 72 hr after CM exposure.

RESULTS: CM volume was higher in the Control group than in the DyeVert group (130 [120–188] mL vs. 99 [69–136] mL; $p < .001$). In the DyeVert group the mean percent CM volume saved was 38 \pm 13%. AKI occurred in 7/90 patients (8%) in the DyeVert group and in 17/90 (19%) patients in the Control group (odds ratio = 0.37; 95% confidence interval 0.14–0.95; $p = .047$).

CONCLUSIONS: This preliminary result suggests that CM volume reduction obtained by the DyeVert™ system is an effective strategy to prevent AKI in ACS patients undergoing invasive procedure.

KEYWORDS
acute coronary syndrome, acute kidney injury, contrast media

1 | INTRODUCTION

Acute kidney injury (AKI) is a common complication in patients suffering from acute coronary syndrome (ACS).^{1,2} This complication has been associated with higher early and late adverse events.^{3,4} The pathogenesis of AKI in ACS patients is multifactorial.^{5,6} Although some data suggest that it is unlikely that contrast media (CM) is a primary pathogenic factor responsible for renal dysfunction in ACS patient treated by percutaneous coronary intervention (PCI), avoiding excess CM volume has been advocated.^{6–8}

The DyeVert™ system (Isprey Medical Inc., Minnesota, MN) is a novel device designed to reduce CM volume during coronary procedures, while maintaining fluoroscopic image quality.^{9,10} The aim of the present study is to assess whether the use of the DyeVert system is


Catheter Cardiovasc Interv. 2020;1–9.

wileyonlinelibrary.com/journal/ccd

© 2020 Wiley Periodicals LLC. | 1

Cost Savings with DyeVert + AKI Reduction

- 4 hospitals
 - DyeVert vs. Control Group Data
- Overall number-needed-to treat
 - 8 to 25 to prevent 1 AKI event
 - Mean NNT = 15
- High cost of an AKI event vs low cost of DyeVert (\$15,000 vs \$350/case)
- Per case cost savings \$250 to \$1,525
- Conservative estimate:
 - Does not consider ongoing renal morbidity
 - Does not consider related downstream health outcomes



Author/Year	Absolute AKI Reduction	Overall NNT (1/AR)	DyeVert Cost (All NNT Cases)	Cost per AKI Event	Cost Offset (Savings Per Case)
Sattar et al, 2018	4.0%	25	\$8,750	\$15,000	\$250
Bunney et al, 2019	5.9%	17	\$5,950	\$15,000	\$532
Briguori et al, 2020	11.0%	9	\$3,150	\$15,000	\$1,317
Kutschman, 2019	12.4%	8	\$2,800	\$15,000	\$1,525
Average	8.3%	15	\$5,163		\$906

- Bunney R, Saenger E, Shah C, Harris S, et al. Contemporary use of contrast dye reduction technology in a tertiary academic hospital: Patient characteristics and acute kidney injury outcomes following percutaneous coronary interventions. Poster presentation at the ACC Quality Summit. 2019.
- Castro D, Dang TT. Reducing Contrast-Induced Acute Kidney Injury in the Cath Lab. Poster presentation at the ACC Quality Summit. 2018. (also published in: Castro D, Marban A. Reducing Contrast-Induced Acute Kidney Injury in the Cath Lab at Houston Methodist Sugar Land Hospital. Cath Lab Digest. 2018;26(6)).
- Kutschman, R. Clinical and Economic Outcomes of a Comprehensive Clinical Quality Initiative for Reducing Acute Kidney Injury in Chronic Kidney Disease Patients Undergoing Coronary Angiography./ JACC. 2019;74(13):Supplement TCT-617.
- Sattar A, Darby M, Schnatz R, El-Hamdani E. Impact of using DyeVert Plus on incidence of acute kidney injury after cardiac catheterizations with coronary interventions in high risk patients. Poster presentation at the West Virginia ACC Chapter meeting. April 2018.
- Thukral N. Kidney Care Protocols in the Cath Lab. TCT 2019 Presentation. Accessed on October 14, 2019 from: <https://www.tctmd.com/slide/aki-reduction-cath-lab-protocol-development-implementation-and-outcomes>.
- Turner, et al. Real-World Impact of a Quality Improvement Program for AKI Prevention in the Cardiac Cath Lab. SCAI Abstract (II-85); Catheter Cardiovasc Interv; Vol 95, Issue S2, May 1, 2020.
- Cameron, et al. Reduction of Contrast Induced Acute Kidney Injury in a Cardiac Catheterization Laboratory: A Quality Improvement Initiative. SCAI Abstract (I-34); Catheter Cardiovasc Interv; Vol 95, Issue S2, May 1, 2020. E-poster abstract presented at the 2020 SCAI Virtual Conference Meeting: <https://virtual2020.scai.org/p/46024>.
- Briguori C, Golino M, Porchetta N, et al. Impact of a contrast media volume control device on acute kidney injury rate in patients with acute coronary syndrome [published online ahead of print, 2020 Jul 18]. *Catheter Cardiovasc Interv*. 2020;10.1002/ccd.29136. doi:10.1002/ccd.29136.

Thank you

*Please visit www.ospreymed.com for
more information*

Indications, contraindications, warnings, and instructions for use can be found in the product labeling supplied with each device. Risks may include air emboli and infection.

*CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a
physician*

be kind to KIDNEYS 