## DyeVert™ Plus EZ Disposable Kit INSTRUCTIONS FOR USE

### DEVICE DESCRIPTION

The Osprey Medical DyeVert™ Plus EZ Contrast Reduction System is a compatible device to manual contrast injections and provides fluid pathway resistance modulation such that excess contrast volume (i.e. contrast that is not needed for diagnostic or therapeutic purposes) is minimized in the patient's vasculature and total contrast agent volume reduction occurs; while maintaining adequate image quality. Age, diabetes, moderate and severe chronic kidney disease and heart failure on presentation are leading factors for when to consider renal protection measures such as contrast minimization tools and processes.

The DyeVert Plus EZ Contrast Reduction System consists of the:

- 1) Display (provided separately) and
- 2) DyeVert Plus EZ Disposable Kit.

The Display may be either the Osprey Medical Contrast Monitoring Display or the Osprey Medical Smart Monitor. The *DyeVert Plus EZ Disposable Kit* consists of the Smart Syringe and DyeVert Plus EZ Module.

The DyeVert Plus EZ Disposable Kit is intended to be used with the Display to allow monitoring and display of contrast volumes manually injected. Volumes are displayed and compared to physician entered contrast usage thresholds during angiographic procedures.

The DyeVert Plus EZ Module has been designed for use with standard injection syringes and manifolds with Luer fittings that have been demonstrated to comply with ISO 594 "Conical fittings with a 6% luer taper for syringes, needles and certain other medical equipment"; and the catheter configurations listed below. Utilization of catheters beyond those listed has not been substantiated.

Diagnostic	Guide	Guide w/Rx	Guide w/OTW
4F	-	-	-
5F	5F	-	-
6F	6F	6F	6F
-	7F	7F	7F

### MODEL NUMBER SELECTION- Disposables

Model Number	Contrast Viscosity Range at 20 C
HV-EZ-XXX-EU	16.3 to 26.6 cps (mPa*s)
HV-EZ-XXX-EU-10	
LV-EZ-XXX-EU	8.8 to 16.3 cps (mPa*s)
LV-EZ-XXX-EU-10	

Note: "XXX" identifier is denoted as available syringe configurations

### CMS CLASSIFICATION

- The DyeVert Plus EZ Disposable Kit is Internally Powered
- Degree of protection against electric shock: Type CF applied part
- Equipment not suitable for use in the presence of flammable mixtures

### APPLIED PARTS

Smart Syringe, DyeVert contrast line and T-Connector of the DyeVert Plus EZ Module and Contrast Collection Bag.

### INTENDED USE

The DyeVert™ Plus EZ Contrast Reduction System is intended to reduce the amount of contrast media administered during procedures requiring the injection of contrast media. Clinical evidence has demonstrated that contrast media can be toxic to the kidneys, leading to Contrast Induced Acute Kidney Injury.

### CONTRAINDICATIONS

Not for use with power injectors.

### WARNINGS

Disposables are for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocess or resterilization may create a risk of patient infection which could lead to injury, illness or death.

The DyeVert Plus EZ Disposable Kit should be connected to single-use contrast source devices only. The DyeVert Plus EZ Disposable Kit does not prevent nor protect against contamination or microbial ingress to or from a single-use contrast source. The DyeVert Plus EZ Disposable Kit and single-use contrast source, including any remaining contrast solution in the single-use contrast source, should be disposed of once the procedure is completed.

Do not use if product packaging appears compromised or damaged.

Do not use catheters or contrast agents not listed in these Instructions for Use (IFU).

Refer to the Display Instructions for Use; and ensure all labeling and Instructions for Use is followed.

- Smart Monitor Instructions for Use and labeling,
- Contrast Monitoring Display Instructions for Use and labeling

Portable and mobile RF communications equipment may affect the devices. The Display should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Display should be observed to verify normal operation.

Use of accessories, transducers and cables other than those specified or provided by Osprey Medical could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Display or disposables, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Please refer to the contrast agent Labeling for dosage recommendations, warnings, contraindications, detail of reported adverse event types and detailed directions for use associated with contrast administration.

For accurate Smart Syringe % contrast concentration values, ensure Smart Syringe is initially primed with 100% contrast media and the contrast source is 100% contrast.

Refer to the Display Instructions for Use for the electromagnetic compliance information of this system.

### PRECAUTIONS

In the event the device malfunctions or changes in performance that is not expected, discontinue use immediately and report experience to Osprey Medical representative. In the event contrast monitoring is not available, the procedure may continue without contrast monitoring and the DyeVert Plus EZ Module will continue to save contrast. If the incident is considered reportable (e.g. serious) by the regulating authority, please ensure the incident is also reported to the qualifying regulatory authority.

Using the DyeVert Plus EZ off-label may result in undesired affects such as poor imaging, lack of contrast reduction or no contrast monitoring.

The DyeVert Plus EZ is designed to be used with non-diluted, room temperature (non-warmed) contrast media only.

As with any device used for injecting contrast media into a patient, care should be taken to ensure all air has been removed from the lines, prior to injection, to avoid air embolization.

Use only light tapping if necessary to remove air while priming the system. Do not use tools (hemostats or other instruments).

During injections, the Smart Syringe, DyeVert Stopcock, and manifold (not provided) should be placed where it is supported on a flat surface across the length of the connections to prevent breakage by unintended bending when manipulating the syringe.

Be cautious to not over-tighten on luer connections when connecting the Smart Syringe or DyeVert Plus EZ Module to the manifold.

The Smart Syringe or DyeVert Plus EZ Module should not be immersed in contrast or saline.

Osprey Medical recommends users follow hospital policy/procedure and physician recommendation on the appropriate total cumulative volume of contrast media used in a patient. The device is not intended to prevent manual injection of contrast media.

User should ensure the Smart Syringe plunger is idle for a minimum of 1/2 second (Dwell Time) prior to switching between contrast & saline when aspirating. If needed, contact Osprey Medical to adjust Smart Syringe plunger dwell time to align with user preference. Inaccurate cumulative volume may be displayed if dwell time is not reached prior to switching between aspirating contrast & aspirating saline.

The graduations on the Contrast Collection Bag are approximate and are not intended for accurate diversion volume measurement. In addition to contrast, the Contrast Collection Bag may contain saline, blood or other fluids.

### POTENTIAL PATIENT ADVERSE EVENTS

Possible adverse effects include but are not limited to: air embolism, infection.

### HOW SUPPLIED

Disposable, sterile system components have been sterilized with ethylene oxide (EO).

EU Declaration of Conformities are available upon request.

#### STORAGE

Store the DyeVert Plus EZ Disposable Kit between -15°C and +38°C (5°F and 100°F).

#### MAINTENANCE and REPAIR

Maintenance is not required.

#### DISPOSAL

Discard Smart Syringe, DyeVert Plus EZ Module with Contrast Collection Bag and single-use contrast source according to hospital procedures. The contrast in the Contrast Collection Bag should not be reused.

The Smart Syringe and DyeVert Plus EZ Module have been designed to comply with battery waste regulations. At the end of device use, the alkaline batteries are inside an infected medical device. In this instance, the Batteries Directive (2006/66/EC), the WEEE Directive and Federal Law DO NOT apply to the collection or waste treatment of the alkaline battery.

Follow local governing ordinances regarding disposal. Do not incinerate as the enclosed batteries may explode at excessive temperatures.

#### OPERATING CONDITIONS

The system is intended to be used in a standard hospital cath lab environment under the following conditions:

Temperature: 10°C to 27°C (50°F to 80°F)

Relative humidity 0% to 85%, noncondensing

The system is not intended be used near active high frequency surgical equipment where the intensity of electromagnetic disturbances is high.

#### MASS of DEVICES

DyeVert Plus EZ Module with Contrast Collection Bag 135g  
Smart Syringe 44g

#### PHYSICIAN TRAINING INFORMATION



Qualified physicians should be knowledgeable with Cath lab procedures, techniques and contrast media usage.

No additional special skills or training are required to operate the system, but physicians should be thoroughly familiar with the DyeVert Plus EZ Contrast Reduction System supporting material including all product labeling. Physicians may contact Osprey Medical to request product review.

#### DIRECTIONS FOR USE


##### System Overview

The DyeVert Plus EZ Disposable Kit is for use with the Monitor through wireless communication. The system allows user input and monitoring of cumulative volume, injection volume and user inputted threshold. A user selectable triple-toggle field allows for display of either the % contrast in Smart Syringe, % Contrast Saved for each injection or % average contrast saved throughout the case. The system with the DyeVert Plus EZ Module allows contrast modulation to the patient resulting in a total contrast volume savings.


#### Modulation System

##### A. Assemble and Prime Manifold

Note: It is recommended to leave all battery pull tabs in the disposable devices until all priming is complete and the case is ready to begin. This will ensure maximum battery life and performance throughout the case.


- 1) Introduce disposables into sterile field per normal procedure.
- 2) Attach contrast collection bag to DyeVert Module Reservoir.
- 3) Attach contrast source line to top of DyeVert Module Reservoir.
- 4) Attach DyeVert contrast line to contrast port on manifold. Ensure swivel luer connection is tightened and securely attached to manifold.
- 5) Attach DyeVert T-Connector to manifold.
- 6) Turn DyeVert 1-Way Stopcock to "OFF" (parallel to manifold) to isolate the DyeVert Module Reservoir.


- 7) Remove air from Smart Syringe and attach to DyeVert T-Connector.
- 8) Prime manifold and all lines per normal procedure.

#### **B. Prime DyeVert Plus EZ Module With Contrast**

- 1) Turn manifold contrast line stopcock to aspiration position.
- 2) Aspirate approximately 6mL of contrast into Smart Syringe.
- 3) Turn manifold contrast line stopcock to "OFF" to Smart Syringe.
- 4) Turn the DyeVert stopcock to "ON" to savings position (parallel to diversion line).

- 5) Inject approximately 4-5 mL of contrast through diversion line until contrast is visible in reservoir. Ensure that all air is clear from diversion line up to the dual check valve. Note: Air may be present beyond the dual check valve.

- 6) Turn contrast stopcock "OFF" to contrast supply

#### **C. Wirelessly Connect Disposables to Display**

Note: Assemble and prime system prior to wirelessly connecting disposables  
 1) **Remove battery pull tab** from the Smart Syringe by pulling the tab straight out from the syringe plunger. **Remove battery pull tab** from the DyeVert Plus EZ Module by pulling the tab straight out from the DyeVert Plus EZ Module (i.e. perpendicular from the back of the Module). Do NOT hold the DyeVert Plus EZ Module down against the table such that the pull tab is against the table when removing the pull tab..

- 2) **Select New Case**. Wirelessly connect Smart Syringe & DyeVert Plus EZ Module to the Display by selecting the device to connect on the Monitor.
- 3) **Confirm wireless connection** of the Display and disposables. Following scanning, the LEDs on the Smart Syringe & DyeVert Plus EZ Module will flash in the same pattern as shown on the Display for proper identification.

- a. Select Yes on the Display or move the plunger of the Smart Syringe > 2 mL to confirm Syringe flashing.

- b. Select Yes on the Display or press Pause Button on DyeVert Plus EZ Module to confirm Module flashing.

Note: Selecting No will result in further scanning.

Note: The green LED on the device will remain on for 10 seconds when the device is connected to the Display.

#### **D. System Operation**

For detailed instructions, refer to the Display Instructions for Use.

Enter the physician specified threshold. If no threshold is desired, enter "0".

**ENSURE ALL PRIMING IS COMPLETED BEFORE PRESSING Start Case.**

- 1) **Select Start Case** to begin contrast accounting.

#### **2) Contrast Injection**

- a. Turn DyeVert Plus EZ Module stopcock ON.

- b. **Aspirate and inject** contrast per normal operating procedure.

**Note:** Per normal injection practice, additional vessel opacity can be achieved by increasing syringe injection speed.

#### **3) Contrast Accounting: Pause and Resume**

The DyeVert Plus EZ Module has a Pause button which allows the user to manually pause/resume contrast accounting.

When the system is in "Pause" mode, cumulative contrast accounting is suspended.

##### **a. To Pause**

###### **Pause Lock Disabled/Enabled**

Select Pause on Display or press DyeVert Plus EZ Module pause button. LED's on Smart Syringe & Module will blink yellow and background on Display will change to yellow.

##### **b. To Resume**

###### **Pause Lock Enabled**

Aspirate contrast into Smart Syringe (resumes automatically), select Resume on Display, or press DyeVert Plus EZ Module pause button. LED's on Smart Syringe and Module will be green for 10 seconds and the Display will return to normal background.

###### **Pause Lock Enabled**

Select Resume on Display, or press DyeVert Plus EZ Module pause button. LED's on Smart Syringe and Module will be green for 10 seconds and the Display will return to normal background.

If contrast is injected back to the source, the system does NOT need to be paused. The system will automatically recognize this as not being administered to the patient. To maintain system accuracy, the user should ensure that only 100% contrast is injected back to the source.

#### **Contrast Accounting Accuracy Tips:**

- Ensure system is **Paused** if contrast is not being injected to patient.
- Ensure system is **Active** (unpaused) when contrast is being injected to patient.
- Ensure the minimum dwell time (1/2 second) is achieved between aspirations of contrast and saline.


#### **4) Contrast Savings Bypass With Contrast Monitoring**

- a. Turn DyeVert Stopcock OFF to Reservoir.
- b. To resume contrast savings turn DyeVert stopcock ON.

#### **5) DRUG ADMINISTRATION**

IF DRUGS ARE TO BE ADMINISTERED TO THE PATIENT THROUGH THE DYEVERT SYSTEM ENSURE THAT THE DYEVERT STOPCOCK IS **OFF** TO SAVINGS. FAILURE TO TURN THE DYEVERT

**STOPCOCK OFF TO SAVINGS COULD RESULT IN A LOWER DOSE OF THE DRUG BEING ADMINISTERED.**

#### DyeVert Plus EZ Contrast Reduction System FCC, IC Information

DyeVert Plus EZ Module FCC ID: 2AHUPDV, IC: 26526-DV  
Smart Syringe FCC ID: 2AHUPSS, IC: 26526-SS

Classifications per IEC 60601-1 / UL 60601-1:

Degree of Protection:

Type CF-Applied Part

Mode of Operation:

Noncontinuous

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

**WARNING** This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions: 1. This device may not cause interference. 2. This device must accept any interference, including interference that may cause undesired operation of the device.

**WARNING:** Changes or modifications to the Display, DyeVert Plus EZ Module or Smart Syringe not expressly approved by Osprey Medical, Inc. could void the user's authority to operate the equipment.

#### ELECTROMAGNETIC INTERFERENCE PRECAUTIONS

This equipment has been tested and found to comply with the limits for a Group 1 Class B device, pursuant to IEC/EN 60601-1-2, 4<sup>th</sup> edition. These limits are designed to provide reasonable protection against harmful interference. This equipment, if not installed and used in accordance with the instructions, may cause harmful interference to other equipment. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the other equipment.
- Increase the separation between the Display and the other equipment.
- Connect the other equipment into an outlet on a circuit different from that to which Display is connected.
- Consult Osprey Medical for help.

**WARNING:** Portable and mobile RF communications equipment may affect the devices. The Display should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Display should be observed to verify normal operation.

**WARNING:** Use of accessories, transducers and cables other than those specified or provided by Osprey Medical could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

**WARNING:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Display or disposables, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

If performance of the system is lost or degraded due to electromagnetic interference, the procedure may continue without contrast monitoring.

The Display and disposables have Bluetooth transceivers using short-wavelength UHF radio waves in the ISM band from 2.4 to 2.485 GHz. The DyeVert Plus EZ Module and Smart Syringe transmit a maximum signal strength of 0.063 mW (-12 dBm ERP) using GSFM modulation as per IEEE 802.15.1 Bluetooth standard and the Bluetooth SIG Working Group specification Version 4.0+.

**WARNING:** Refer to the Smart Monitor or Contrast Monitoring Display IFU for the electromagnetic compliance information of this system.

**DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY**

THERE IS NO WARRANTY EXPRESSED OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR OF FITNESS FOR A PARTICULAR PURPOSE, ON THE OSPREY MEDICAL PRODUCT(S) DESCRIBED IN THIS PUBLICATION. IN THE EVENT OF ANY DEFECT OR NONCONFORMITY OF OR TO THIS PRODUCT(S), OSPREY MEDICAL'S LIABILITY SHALL NOT BE IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT(S) TO BUYER. UNDER NO CIRCUMSTANCES SHALL OSPREY MEDICAL BE LIABLE FOR ANY DIRECT, INDIRECT OR CONSEQUENTIAL DAMAGES BASED UPON BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, STRICT TORT OR ANY OTHER THEORY ARISING OUT OF THE PURCHASE, USE OR REUSE OF THIS PRODUCT(S). OSPREY MEDICAL NEITHER ASSUMES NOR AUTHORIZES ANY PERSON TO ASSUME FOR IT ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH OSPREY MEDICAL PRODUCT(S). Descriptions or specifications in Osprey Medical printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Packaging Symbol Definitions

	Expiration Date YYYY-MM. Use by last day of month (MM).		Manufacturer	<b>STERILE   EO</b>	Sterilized by ethylene oxide
	Consult electronic Instructions for Use		Keep Dry	<b>REF</b>	Model number
	Single Use		European Conformity	<b>LOT</b>	Lot Number
	Prescription Only		Do Not Use if Package is Damaged		EN IEC 60601-1 applied part Type CF Defibrillation-proof
	Medical Device		Consult instructions for use for important information		MR Unsafe
	Temperature Limit		Identification of Sterile Barrier		Refer to
	DyeVert Plus EZ Module		Syringe		
	EU authorized representative		Australian Communications and Media Authority mark		

Osprey Medical Inc.  
5600 Rowland Road, Suite 250  
Minnetonka, MN 55343  
USA  
Customer Service Toll-Free:  
1-855-860-7584 Fax: 1-855-883-4365  
customerservice@ospreymed.com  
www.ospreymed.com

Osprey & DyeVert are trademarks of Osprey Medical Inc.  
©Osprey Medical, Inc. 2020. All Rights Reserved.

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Australian Sponsor  
Osprey Medical, Pty  
Level 13, 41 Exhibition Street  
Melbourne, Victoria 3000 Australia

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MedPass SAS