## DyeTect™ Contrast Monitoring Disposable Kit INSTRUCTIONS FOR USE

**CAUTION: Federal (USA) Law restricts this device to sale by or on the order of a physician.**

### DEVICE DESCRIPTION

The DyeTect™ Contrast Monitoring Disposable Kit consists of a Smart Syringe and Pressure Module to be used with the Display during controlled infusion for procedures requiring injection of contrast media. The DyeTect Contrast Monitoring System allows for monitoring and display of contrast volumes manually injected.

The Pressure Module has been designed for use with manifolds with Luer fittings that have been demonstrated to comply with ISO 594 "Conical fittings with a 6% luer taper for syringes, needles and certain other medical equipment".

### MODEL NUMBER SELECTION - Disposables

Model Number	Syringe Description
RRL	Locking Luer, Ring Plunger – Ring Grip
PRL	Locking Luer, Palm Plunger – Ring Grip
PWL	Locking Luer, Palm Plunger – Wing Grip
RRS	Swivel Luer, Ring Plunger – Ring Grip
PRS	Swivel Luer, Palm Plunger – Ring Grip

### CMS CLASSIFICATION

- Type of protection against electric shock: Class 1.
- Degree of protection against electric shock: Type CF applied part
- Equipment not suitable for use in the presence of flammable mixtures

### APPLIED PARTS

Smart Syringe and Pressure Transducer of the DyeTect Pressure Module. All parts are single use, disposables.

### INTENDED USE

The Contrast Monitoring Disposable Kit consists of a Smart Syringe and a Pressure Module to be used with the Display during controlled infusion for procedures requiring injection of contrast media. The Osprey Medical DyeTect Contrast Monitoring System allows for real-time monitoring and display of contrast volumes manually injected.

### INDICATION FOR USE

The DyeTect™ Contrast Monitoring Disposable Kit consists of a Smart Syringe and Pressure Module to be used with the Display during angiographic or CT procedures requiring controlled infusion of radiopaque contrast media.

### CONTRAINdicATIONS

Not for use with power injectors.

### WARNINGS

Disposables are for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocess or resterilization may create a risk of patient infection which could lead to injury, illness or death.

The DyeTect Contrast Monitoring Disposable Kit should be connected to single-use contrast source devices only. The Pressure Module does not prevent nor protect against contamination or microbial ingress to or from a single-use contrast source. The DyeTect Contrast Monitoring Disposable Kit and single-use contrast source, including any remaining contrast solution in the single-use contrast source, should be disposed of following the procedure.

Do not use if product packaging appears compromised.

Refer to the specific Display Instructions for Use and ensure all labeling and Instructions for Use are followed:

- Smart Monitor Instructions for Use and labeling, or
- Contrast Monitoring Display Instructions for Use and labeling

Please refer to the contrast agent Labeling for dosage recommendations, warnings, contraindications, detail of reported adverse event types and detailed directions for use associated with contrast administration.

For accurate Smart Syringe % contrast concentration values, ensure Smart Syringe is initially primed with 100% contrast media and the contrast source is 100% contrast.

Changes or modifications to the Monitor, Smart Syringe or Pressure Module not expressly approved by Osprey Medical could void the user's authority to operate the equipment.

Portable and mobile RF communications equipment may affect the devices. The Monitor should not be used adjacent to or stacked with other equipment. If

adjacent or stacked use is necessary, the Display should be observed to verify normal operation.

Use of accessories, transducers and cables other than those specified or provided by Osprey Medical could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Monitor or disposables, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Refer to the Smart Monitor or Contrast Monitoring Display Instructions for Use for the electromagnetic compliance information of this system.

### PRECAUTIONS

In the event the device malfunctions or changes in performance that is not expected, discontinue use immediately and report experience to Osprey Medical representative. In the event contrast monitoring is not available, the procedure may continue without contrast monitoring. If the incident is considered reportable (e.g. serious) by the regulating authority, please ensure the incident is also reported to the qualifying regulatory authority.

Using the DyeTect off-label may result in undesired affects such as no contrast monitoring.

As with any device used for injecting contrast media into a patient, care should be taken to assure all air has been removed from the lines, prior to injection, to avoid air embolization.

Use only light tapping if necessary to remove air while priming the system. Do not use tools (hemostats or other instruments).

Be cautious to not over-tighten on luer connections when connecting the Smart Syringe and/or Pressure Module to a manifold.

The Smart Syringe or Pressure Module should not be immersed in contrast or saline.

Osprey Medical recommends users follow hospital policy/procedure and physician recommendation on the appropriate total cumulative volume of contrast media used in a patient. The device is not intended to prevent manual injection of contrast media.

User should ensure the Smart Syringe plunger is idle for a minimum of 1/2 second (Dwell Time) prior to switching between contrast & saline when aspirating. If needed, contact Osprey Medical to adjust Smart Syringe plunger dwell time to align with user preference. Inaccurate cumulative volume may be displayed if dwell time is not reached prior to switching between aspirating contrast & aspirating saline.

### POTENTIAL PATIENT ADVERSE EVENTS

Possible adverse effects include but are not limited to: air embolism, infection.

### HOW SUPPLIED

Disposable, sterile system components have been sterilized with ethylene oxide (EO).

EU Declaration of Conformities are available upon request.

### STORAGE

Store the DyeTect Disposable Kit between -15°C and +38°C (5°F and 100°F).

### MAINTENANCE and REPAIR

Maintenance is not required.

### DISPOSAL

Discard Smart Syringe, Pressure Module and single-use contrast source according to hospital procedures.

The Smart Syringe and Pressure Module have been designed to comply with battery waste regulations. At the end of device use, the alkaline batteries are inside an infected medical device. In this instance, the Batteries Directive (2006/66/EC), the WEEE Directive and Federal Law DO NOT apply to the collection or waste treatment of the alkaline battery. Follow local governing ordinances regarding disposal. Do not incinerate as the enclosed batteries may explode at excessive temperatures.

### OPERATING CONDITIONS

The system is intended to be used in a standard hospital Cath lab environment under the following conditions:

Temperature: 10°C to 27°C (50°F to 80°F)

Relative humidity 0% to 85%, noncondensing

The system is not intended to be used near active high frequency surgical equipment where the intensity of electromagnetic disturbances is high.

**MRI SAFTEY INFORMATION**

The Smart Syringe and DyeTect Pressure Module are MR UNSAFE. Keep them outside the MRI Scanner room.

**MASS of DEVICES**

Smart Syringe	44g
Pressure Module	38g

**PHYSICIAN TRAINING INFORMATION**

Qualified physicians should be knowledgeable with Cath lab procedures, techniques and contrast media usage.


No additional special skills or training are required to operate the system but physicians should be thoroughly familiar with the DyeTect Contrast Monitoring System supporting material including all product labeling. Physicians may contact Osprey Medical to request product review.

**DIRECTIONS FOR USE**

Note: For Detailed instructions for the Display, refer to the Contrast Monitoring Display Instructions for Use or the Smart Monitor Instructions for Use as applicable.

**System Overview**

The DyeTect Contrast Monitoring Disposable Kit is for use with the Display through wireless communication. The system allows user input and monitoring of cumulative volume, injection volume, % contrast concentration in Smart Syringe and user inputted threshold throughout the case.

- 3) **Confirm wireless connection** of the Display and disposables. Following scanning, the LED's on the Smart Syringe & Pressure Module will flash in the same pattern as shown on the Display for proper identification.
  - a) Select Yes on the Display or move the plunger of the Smart Syringe > 2 mL to confirm Syringe flashing.
  - b) Select Yes on the Display or press Pause Button on Pressure Module to confirm Module flashing.

Note: Selecting No will result in further scanning and will not allow connection of the rejected device to the Monitor.

Note: The green LED on the syringe and Module, respectively, will remain on for 10 seconds when each device is connected to the Display.

#### Start Case

- 1) Enter the physician specified threshold. If no threshold is desired, enter "0".

#### **ENSURE ALL PRIMING IS COMPLETED BEFORE PRESSING Start Case ON THE Display .**

For Detailed instructions, refer to the Contrast Monitoring Display or Smart Monitor Instructions for Use.

- 2) **Select Start Case** on Display to begin contrast accounting.
- 3) **Perform procedure**

#### Pause/Resume

The DyeTect Pressure Module has a Pause button which allows the user to manually pause/resume contrast accounting.

When the system is in "Pause" mode, cumulative contrast accounting is suspended. "Resume" mode allows contrast accounting

- 1) **To Pause** select **Pause** on Display or press DyeTect Module pause button. LED's on Smart Syringe & Module will blink yellow and background on Display will change to yellow.

- 2) **To Resume**

Pause Lock Disabled

Aspirate contrast into Smart Syringe (resumes automatically) or select **Resume** on Display or press Module pause button. LED's on Smart Syringe and Module will be green for 10 seconds and the Monitor will return to normal background.

Pause Lock Enabled

Select **Resume** on Display or press Module pause button. LED's on Smart Syringe and Module will be green for 10 seconds and the Monitor will return to normal background.

If contrast is injected back to the source, the system does NOT need to be paused. The system will automatically recognize this as not being administered to the patient. To maintain system accuracy, the user should ensure that only 100% contrast is injected back to the source.

#### Contrast Accounting Accuracy Tips:

- Ensure system is Paused if contrast is not being injected to patient.
- Ensure system is Active (unpaused) when contrast is being injected to patient.
- Ensure the minimum dwell time (1/2 second) is achieved between aspirations of contrast and saline.

#### End Case

- 1) **Press End Case button**
- 2) When prompted, confirm that you want to end the case

**Caution:** Ending case will permanently disable communication between the Monitor and disposables.


Case summary screen appears showing cumulative volume of contrast administered to patient and the percentage of physician specified threshold.

The Smart Syringe and Pressure Module contain batteries. Discard the Smart Syringe and Pressure Module and single-use contrast source according to hospital procedures.

Follow local governing ordinances regarding disposal. Do not incinerate as the enclosed batteries may explode at excessive temperatures

Note: To view summaries for previous cases, from **Main Menu** select **Settings, User Settings, Case History**.

**Select Main Menu** to return to the main menu and start a new case.

The Pressure Module connects wirelessly to the Monitor.

It includes a module with a pause button and a pressure transducer.

The Smart Syringe and Pressure Module have two LEDs to provide user feedback.

- **Solid green** (active for 10 seconds) indicates injections will be counted towards cumulative volume
- **Flashing Yellow** indicates system is paused and injections will not be counted towards cumulative volume

#### Assemble and Prime

Note: It is recommended to leave all battery pull tabs in the disposable devices until all priming is complete and the case is ready to begin. This will ensure maximum battery life and performance throughout the case.

- 1) Introduce disposables into sterile field per normal procedure.
- 2) Attach contrast source line to the female luer located in the pressure transducer of the Pressure Module.
- 3) Attach the male luer of the pressure transducer to the Manifold, and the Manifold to the Smart Syringe.
- 4) Prime manifold and all lines per normal procedure.

#### Wirelessly Connect Disposables to Monitor

Note: Assemble and prime system prior to wirelessly connecting disposables

- 1) **Remove battery pull tab** battery pull tab from the Smart Syringe by pulling the tab straight out from the syringe plunger. **Remove battery pull tab** from the Pressure Module by pulling the tab straight out from the Module (i.e. perpendicular from the back of the Module). Do NOT hold the Pressure Module down against the table such that the pull tab is against the table when removing the pull tab.
- 2) **Select New Case.** Wirelessly connect Smart Syringe & Pressure Module to the Monitor by selecting the device to connect on the display.

### System Discontinuation

**Smart Monitor:** Power Down by pressing and holding the power button on the end of the Monitor, then swiping the screen as indicated.

**Contrast Monitoring Display:** Select Power Down from the main menu and then select YES.

### Contrast Monitoring System FCC, IC Information

Smart Syringe FCC ID: 2AHUPSS, IC: 26526-SS  
DyeTect Module FCC ID: 2AHUPPM, IC: 26526-PM

Classifications per IEC 60601-1 / UL 60601-1:

Degree of Protection:	Type CF-Applied Part
Mode of Operation:	Noncontinuous

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

**WARNING** This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions: 1. This device may not cause interference. 2. This device must accept any interference, including interference that may cause undesired operation of the device.

**WARNING:** Changes or modifications to the Monitor, Smart Syringe or Pressure Module not expressly approved by Osprey Medical could void the user's authority to operate the equipment.

### ELECTROMAGNETIC INTERFERENCE PRECAUTIONS

This equipment has been tested and found to comply with the limits for a Group 1 Class B device, pursuant to IEC/EN 60601-1-2, 4<sup>th</sup> edition. These limits are designed to provide reasonable protection against harmful interference. This equipment, if not installed and used in accordance with the instructions, may cause harmful interference to other equipment. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the other equipment.
- Increase the separation between the Monitor and the other equipment.
- Connect the other equipment into an outlet on a circuit different from that to which display is connected.
- Consult Osprey Medical for help.

**WARNING:** Portable and mobile RF communications equipment may affect the devices. The Monitor should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Display should be observed to verify normal operation.

**WARNING:** Use of accessories, transducers and cables other than those specified or provided by Osprey Medical could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

**WARNING:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Monitor or disposables, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

If performance of the system is lost or degraded due to electromagnetic interference, the procedure may continue without contrast monitoring.

The Monitor and disposables have Bluetooth transceivers using short-wavelength UHF radio waves in the ISM band from 2.4 to 2.485 GHz. The Pressure Module and Smart Syringe transmit a maximum signal strength of 0.063 mW (-12 dBm ERP) using GSFM modulation as per IEEE 802.15.1 Bluetooth standard and the Bluetooth SIG Working Group specification Version 4.0+.

**WARNING:** Refer to the Smart Monitor or Contrast Monitoring Display Instructions for Use for the electromagnetic compliance information of this system.

### DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

THERE IS NO WARRANTY EXPRESSED OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR OF FITNESS FOR A PARTICULAR PURPOSE, ON THE OSPREY MEDICAL PRODUCT(S) DESCRIBED IN THIS PUBLICATION. IN THE EVENT OF ANY DEFECT OR NONCONFORMITY OF OR TO THIS PRODUCT(S), OSPREY MEDICAL'S LIABILITY SHALL NOT BE IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT(S) TO BUYER. UNDER NO CIRCUMSTANCES SHALL OSPREY MEDICAL BE LIABLE FOR ANY DIRECT, INDIRECT OR CONSEQUENTIAL DAMAGES BASED UPON BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, STRICT TORT OR ANY OTHER THEORY ARISING OUT OF THE PURCHASE, USE OR REUSE OF THIS PRODUCT(S). OSPREY MEDICAL NEITHER ASSUMES NOR AUTHORIZES ANY PERSON TO ASSUME FOR IT ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH OSPREY MEDICAL PRODUCT(S). Descriptions or specifications in Osprey Medical printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Packaging Symbol Definitions					
	Expiration Date YYYY-MM-DD. Use by last day of month (MM).		Manufacturer		Sterilized by ethylene oxide
	Consult electronic Instructions for Use		Keep Dry		Model number
	Single Use		European Conformity		Lot Number
	Prescription Only		Do Not Use if Package is Damaged		EN IEC 60601-1 applied part Type CF Defibrillation- proof
	MR Unsafe		Medical Device		Consult instructions for use for important information
	Identification of Sterile Barrier		Temperature Limit		Syringe
	DyeTect Module		Refer to		
	EU authorized representative		Australian Communications and Media Authority mark		

Osprey Medical, Inc.  
5600 Rowland Road, Suite 250  
Minnetonka, MN 55343  
USA

Customer Service Toll-Free:  
1-855-860-7584 Fax: 1-855-883-4365  
customerservice@ospreymed.com  
www.ospreymed.com

Icon Limited  
South County Business Park, Leopardstown  
Dublin 18, D18 X5R3, Ireland  
SRN IE-AR-000006507

Australian Sponsor  
Osprey Medical, Pty  
Level 13, 41 Exhibition Street  
Melbourne, Victoria 3000 Australia

Osprey & DyeTect are trademarks of Osprey Medical Inc.  
©Osprey Medical Inc. 2020. All Rights Reserved.

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Copyright © Osprey Medical Inc. 2020.