## Smart Syringe INSTRUCTIONS FOR USE

**CAUTION: Federal (USA) Law restricts this device to sale by or on the order of a physician.**

### DEVICE DESCRIPTION

The Smart Syringe is a component to the DyeVert™ Plus Contrast Reduction System (The DyeVert™ Plus Contrast Reduction System has not been licensed in accordance with Canadian law), DyeVert Plus™ EZ Contrast Reduction System and DyeTect™ Contrast Monitoring System (The DyeTect™ Contrast Monitoring System has not been licensed in accordance with Canadian law). The Smart Syringe is not intended to be used independently. The Smart Syringe is a control syringe with wireless communication capability to the Osprey Medical Display.

### MODEL NUMBER SELECTION

Model Number	Syringe Description
R-RRL	Locking Luer, Ring Plunger – Ring Grip
R-PRL	Locking Luer, Palm Plunger – Ring Grip
R-PWL	Locking Luer, Palm Plunger – Wing Grip
R-RRS	Swivel Luer, Ring Plunger – Ring Grip
R-PRS	Swivel Luer, Palm Plunger – Ring Grip

### CMS CLASSIFICATION

- Type of protection against electric shock: Class 1.
- Degree of protection against electric shock: Type CF applied part
- Equipment not suitable for use in the presence of flammable mixtures

### APPLIED PARTS

Smart Syringe

### INTENDED USE

The Smart Syringe is intended to be used for the controlled infusion of radiopaque contrast media for angiographic procedures with the DyeVert Plus and DyeVert Plus EZ Contrast Reduction Systems or the DyeTect Contrast Monitoring System.

### INDICATION FOR USE

The Smart Syringe is to be used with the DyeVert Plus and DyeVert Plus EZ Contrast Reduction Systems or the DyeTect Contrast Monitoring System. The DyeVert Plus and DyeVert Plus EZ Contrast Reduction Systems are to be used for contrast volume reduction and for the monitoring of radiopaque contrast media during angiographic or CT procedures with the following agents: Iodixanol 270 or 320 mgI/mL, Iohexol 300 or 350 mgI/mL and Iopamidol 370 mgI/mL. The DyeTect Contrast Monitoring System is to be used during angiographic or CT procedures requiring controlled infusion of radiopaque contrast media.

### CONTRAINdicATIONS

Not for use with power injectors.

### WARNINGS

Refer to the DyeVert Plus, DyeVert Plus EZ Contrast Reduction Systems or the DyeTect Contrast Monitoring System Instructions for Use for additional information, prior to use.

If the Smart Syringe is not paired with the Contrast Monitoring Display in conjunction with a DyeVert Plus Module, DyeVert Plus EZ Module or DyeTect Pressure Module, accurate accounting of contrast use will not occur.

The Smart Syringe is for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocess or resterilization may create a risk of patient infection which could lead to injury, illness or death.

Do not use if product packaging appears compromised or damaged.

Please refer to the contrast agent Instructions for Use and Labeling for dosage recommendations, warnings, contraindications, detail of reported adverse event types and detailed directions for use associated with contrast administration.

For accurate Smart Syringe % contrast concentration values, ensure Smart Syringe is initially primed with 100% contrast media and the contrast source is 100% contrast.

Portable and mobile RF communications equipment may affect the devices. The Display should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Display should be observed to verify normal operation.

Use of accessories, transducers and cables other than those specified or provided by Osprey Medical could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Display or disposables, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Refer to the Smart Monitor or Contrast Monitoring Display Instructions for Use for the electromagnetic compliance information of this system.

In the event contrast monitoring is not available, the procedure may continue without contrast monitoring.

### PRECAUTIONS

In the event the device malfunctions or changes in performance that is not expected, discontinue use immediately and report experience to Osprey Medical representative. If the incident is considered reportable (e.g. serious) by the regulating authority, please ensure the incident is also reported to the qualifying regulatory authority.

Using the Smart Syringe off-label may result in undesired affects such as poor imaging, lack of contrast reduction or no contrast monitoring.

As with any device used for injecting contrast media into a patient, care should be taken to assure all air has been removed, prior to injection, to avoid air embolization.

Use only light tapping if necessary to remove air while priming the Smart Syringe. Do not use tools (hemostats or other instruments).

Be cautious to not over-tighten on luer connection when connecting the Smart Syringe to a manifold.

The Smart Syringe should not be immersed in contrast or saline.

Osprey Medical recommends users follow hospital policy/procedure and physician recommendation on the appropriate total cumulative volume of contrast media used in a patient. The Smart Syringe is not intended to prevent manual injection of contrast media.

User should ensure the Display dwell time is reached prior to Smart Syringe aspiration and if needed, contact Osprey Medical to adjust dwell time to align with user preference. Inaccurate cumulative volume may be displayed if aspiration occurs prior to dwell time is reached.

### POTENTIAL PATIENT ADVERSE EVENTS

Possible adverse effects include but are not limited to: air embolism, infection.

### HOW SUPPLIED

Disposable, sterile system components have been sterilized with ethylene oxide (EO).

### STORAGE

Store the Smart Syringe between -15°C and +38°C (5°F and 100°F).

### MAINTENANCE and REPAIR

Maintenance is not required.

### DISPOSAL

Discard Smart Syringe according to hospital procedures.

The Smart Syringe has been designed to comply with battery waste regulations. At the end of device use, the alkaline batteries are inside an infected medical device. In this instance, the Batteries Directive (2006/66/EC), the WEEE Directive and Federal Law DO NOT apply to the collection or waste treatment of the alkaline battery.

Follow local governing ordinances regarding disposal. Do not incinerate as the enclosed batteries may explode at excessive temperatures.

### OPERATING CONDITIONS

The system is intended to be used in a standard hospital Cath lab environment under the following conditions:

Temperature: 10°C to 27°C (50°F to 80°F)

Relative humidity 0% to 85%, noncondensing

Pressure = 60 kPa to 106 kPa

The system is not intended to be used near active high frequency surgical equipment where the intensity of electromagnetic disturbances is high.

### MRI SAFETY INFORMATION

The Smart Syringe is MR UNSAFE. Keep it outside the MRI Scanner room.

### MASS of DEVICES

Smart Syringe

44g

### PHYSICIAN TRAINING INFORMATION

Qualified physicians should be knowledgeable with Cath lab procedures, techniques and contrast media usage.

## DIRECTIONS FOR USE

The Smart Syringe is to be used with the DyeVert Plus, DyeVert Plus EZ Contrast Reduction Systems or the DyeTect Contrast Monitoring System. This includes the Contrast Monitoring Display or Smart Monitor and a disposable Module (DyeVert Plus Module, DyeVert Plus EZ Module or Pressure Module). Refer to the DyeVert Plus Contrast Reduction System, DyeVert Plus EZ Contrast Reduction System or the DyeTect Contrast Monitoring System Instructions for Use for additional information, prior to use.

The Smart Syringe has two LEDs to provide user feedback.

1. **Solid green** (active for 10 seconds) indicates injections will be counted towards cumulative volume
2. **Flashing Yellow** indicates system is paused and injections will not be counted towards cumulative volume

### Wirelessly Connect the Smart Syringe to the Contrast Monitoring Display

1. Remove orange battery pull tab from the Smart Syringe.
2. Select **New Case**. Following scanning, the LEDs on the Smart Syringe will flash in the same pattern as shown on the Display for proper identification.

#### Confirm wireless connection

- a) **Select Yes or move the plunger of the Smart Syringe > 2 mL** to confirm Syringe flashing.

Note: Selecting **No** will result in further scanning.

Note: The green LED on the device will remain on for 10 seconds once the device is connected to the Display.

**ENSURE ALL PRIMING IS COMPLETED BEFORE PRESSING Start Case.**

**ENSURE** to wirelessly connect a DyeVert Plus Module or a Pressure Module to the Display.

- b) **Select Start Case** to begin contrast accounting.
- c) **Perform procedure**

**Caution:** Ending case will permanently disable communication between the Display and Smart Syringe.

#### Smart Syringe FCC, IC Information

Smart Syringe FCC, IC: 2AHUPSS, IC: 26526-SS

Classifications per IEC 60601-1 / UL 60601-1:

Degree of Protection:	Type CF-Applied Part
Mode of Operation:	Noncontinuous

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

**WARNING** This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions: 1. This device may not cause interference. 2. This device must accept any interference, including interference that may cause undesired operation of the device.

**WARNING:** Changes or modifications to the Smart Syringe not expressly approved by Osprey Medical could void the user's authority to operate the equipment.

#### SYSTEM ELECTROMAGNETIC INTERFERENCE PRECAUTIONS

This equipment has been tested and found to comply with the limits for a Group 1 Class B device, pursuant to IEC/EN 60601-1-2, 4<sup>th</sup> edition. These limits are designed to provide reasonable protection against harmful interference. This equipment, if not installed and used in accordance with the instructions, may cause harmful interference to other equipment. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

1. Reorient or relocate the other equipment.
2. Increase the separation between the display and the other equipment.
3. Connect the other equipment into an outlet on a circuit different from that to which display is connected.
4. Consult Osprey Medical for help.

**WARNING:** Portable and mobile RF communications equipment may affect the devices. The Display should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Display should be observed to verify normal operation.

**WARNING:** Use of accessories, transducers and cables other than those specified or provided by Osprey Medical could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

**WARNING:** Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Display or disposables, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

If performance of the system is lost or degraded due to electromagnetic interference, the procedure may continue without contrast monitoring.

The Display and disposables have Bluetooth transceivers using short-wavelength UHF radio waves in the ISM band from 2.4 to 2.485 GHz. The Pressure Module and Smart Syringe transmit a maximum signal strength of 0.063 mW (-12 dBm ERP) using GSFM modulation as per IEEE 802.15.1 Bluetooth standard and the Bluetooth SIG Working Group specification Version 4.0+.

**WARNING:** Refer to the Smart Monitor or Contrast Monitoring Display Instructions for Use for the electromagnetic compliance information of this system.

#### DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

THERE IS NO WARRANTY EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR OF FITNESS FOR A PARTICULAR PURPOSE, ON THE OSPREY MEDICAL PRODUCT(S) DESCRIBED IN THIS PUBLICATION. IN THE EVENT OF ANY DEFECT OR NONCONFORMITY OF OR TO THIS PRODUCT(S), OSPREY MEDICAL'S LIABILITY SHALL NOT BE IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT(S) TO BUYER. UNDER NO CIRCUMSTANCES SHALL OSPREY MEDICAL BE LIABLE FOR ANY DIRECT, INDIRECT OR CONSEQUENTIAL DAMAGES BASED UPON BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, STRICT TORT OR ANY OTHER THEORY ARISING OUT OF THE PURCHASE, USE OR REUSE OF THIS PRODUCT(S). OSPREY MEDICAL NEITHER ASSUMES NOR AUTHORIZES ANY PERSON TO ASSUME FOR IT ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH OSPREY MEDICAL PRODUCT(S). Descriptions or specifications in Osprey Medical printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Packaging Symbol Definitions				
	Expiration Date YYYY-MM. Use by last day of month (MM).		Manufacturer	Sterilized by ethylene oxide
	Consult electronic Instructions for Use		Keep Dry	Model number
	Single Use		EN IEC 60601-1 applied part Type CF Defibrillation-proof	Lot Number
	Prescription Only		Do Not Use if Package is Damaged	MR Unsafe
	Medical Device		Consult instructions for use for important information	Identification of Sterile Barrier
	Temperature Limit		Syringe	Refer to
	European Conformity		EU Authorized Representative	

Osprey Medical Inc.

5600 Rowland Road, Suite 250  
Minnetonka, MN 55343

USA

Customer Service Toll-Free:

1-855-860-7584 Fax: 1-855-883-4365

[customerservice@ospreymed.com](mailto:customerservice@ospreymed.com)

[www.ospreymed.com](http://www.ospreymed.com)

Icon Limited

South County Business Park, Leopardstown  
Dublin 18, D18 X5R3, Ireland  
SRN IE-AR-000006507

Australian Sponsor

Osprey Medical, Pty  
Level 13, 41 Exhibition Street  
Melbourne, Victoria 3000 Australia

Osprey, DyeVert and DyeTect are trademarks of Osprey Medical Inc.  
©Osprey Medical Inc. 2020. All Rights Reserved.

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